

NEW PRODUCTS

Digital tape interface

R.C. Electronics (Santa Barbara, CA) introduces the DTI-512 digital tape interface for the Metrum RSR 512 HSDO port to provide high-speed transfer of data to a PC. It has a 1:1 transfer rate to a PC RAM disk, a 4:1 rate to PC hard disk, and a 2:1 rate to HP9000 series 700 workstations. The DTI-512 software screens simulate the RSR 512 recorder format with function keys.


Write in 980 on Reader Service Card.



Gel dryer

The Scienceware gel dryer from Bel-Art Products (Pequannock, NJ) produces secure polyacrylamide and other types of gels in any size up to 11 x 6.5 inches. This device can be used with a heat source, such as a hot plate, for quick drying or without a heat source for slow drying. The unit allows for vacuum sublimation of moisture from the gel if a vacuum pump is attached.

Write in 981 on Reader Service Card.



Instant camera

The new MicroCam SLR instant camera from Polaroid (Montigny-le-Bretonneux, France) is designed for light microscopists. It features a multi-element 10X magnification lens, microprocessor-controlled exposure and filtration, digital light measurement, and a rotary shutter to create instant color and black-and-white prints.

Write in 972 on Reader Service Card.



Diluter

The ASAP gravimetric diluter from Don Whitely Scientific (Shipley, U.K.) automatically prepares samples for plating, ensuring accurate dilution. The user has only to specify the dilution factor and place the sample on the integral balance where it is weighed for the proper automatic dilution.

Write in 974 on Reader Service Card.



Disc dispenser

Unipath (Basingstoke, U.K.) introduces a new ergonomically designed disc dispenser to complement the Oxoid antimicrobial susceptibility testing system. Self-tamping, one-handed operation makes for easy use. The dispenser is available in a 6 or 8 disc format.

Write in 975 on Reader Service Card.



Electrophoresis kit

The new eCAP SDS 14-200 kit from Beckman Instruments (Fullerton, CA) automates SDS protein purity determinations on the P/ACE capillary electrophoresis system. The kit includes a coated capillary, gel buffer, test mix, internal standard-orange G dye, and sample buffer.

Write in 983 on Reader Service Card.



**INTERNATIONAL ADVERTISING
MANAGER (NEW YORK)**
Stephanie J. Nolan 1 (800) 989-1112

**INTERNATIONAL ADVERTISING
MANAGER (LONDON)**
Angela Kays (71) 872 0103

EAST COAST & MIDWEST
Richard T. Kobel 1 (212) 477-3626
65 Bleecker Street
New York, NY 10012
1 (800) 989-1112
FAX: 1 (212) 505-1364

WEST COAST & SOUTHWEST
Lynne Stickrod 1 (415) 781-6422
582 Market Street, Ste. 1407
San Francisco, CA 94104
1 (800) 229-9758
FAX: 1 (415) 781-3805

CANADA
Peter Drake 1 (416) 690-2423
17 Pine Crescent
Toronto, Ontario M4E 1L1
FAX: 1 (416) 868-1213

UNITED KINGDOM
Angela Kays (71) 872-0103
Kathryn Wayman
4 Little Essex Street
London, WC2R 3LF
FAX: (71) 240-2408

GERMANY & AUSTRIA
Sabine Furst (89) 52 40 81
Sandstrasse 41
D80335 München
FAX: (89) 523 22 22

SCANDINAVIA
Andrew Karnig (8) 6 44 00 05
Andrew Karnig & Associates
Finnbodavagen 1
13131 NACKA Sweden
FAX: (8) 6 42 31 50

**FRANCE & HOLLAND
& BELGIUM**
Patricia Hummel (1) 43-87-42-17
5 Rue Joseph Sansboeuf
75008 Paris France
FAX: (1) 43-87-42-15

SWITZERLAND
Verena Loewenthal-Hugi (1) 750 55 11
Interpress Postfach
CH-8103 Unterengstringen
FAX: (1) 750 22 10

JAPAN
Phillip Hamill (3) 267-8751
Shin-Mitsuke Building (4F)
3-6 Ichigaya Tamachi
Shinjuku-ku Tokyo 162
TEL: (3) 267 8751
FAX: (3) 267- 8746

AUSTRALIA
Anne Murrey (2) 299-5677
World Media Network Pty, Level 1
89 York Street, Sydney NSW 2000
FAX: (2) 299-6178

ISRAEL
Heiko Sieger, Intermedia
76 King George St.
PO Box 11820
Tel Aviv, Israel 61116
TEL: (3) 5287815
FAX: (3) 299259